

ENDOPATH

All the proven benefits of ENDOPATH XCEL[®] Trocars with **smudge mitigation**

ENDOPATH XCEL[®] Trocars with OPTIVIEW[®] Technology

Designed for camera ports, OPTIVIEW Technology uses wiping, wicking and absorbing to reduce smudging during endoscope reinsertions

- Reduces rewiping of the endoscope¹
- Minimizes the cumulative time needed to clean the endoscope¹
- Reduces smudging upon endoscope reinsertions more effectively¹
 - At least **59% more effective** than Covidien Versaport[™] and VersaStep[™] Bladeless 5mm Trocars
 - At least **58% more effective** than Applied Medical Kii[®] and Kii Fios[®] 5mm Trocars
 - At least **66% more effective** than Covidien Versaport and VersaStep Bladeless 12mm Trocars
 - At least **74% more effective** than Applied Medical Kii and Kii Fios 12mm Trocars

Low insertion force

Provides controlled and predictable entry

Abdominal wall retention

Integrated thread design on Bladeless, Universal Sleeve and Dilating Tip cannulas help minimize trocar slip-outs

Low drag force

Housing and seal are designed to reduce instrument hang-ups and provide low instrument insertion/extraction forces

Direct visualization

Bladeless tip trocar allows for endoscopic optical entry

ENDOPATH XCEL Trocar with OPTIVIEW Technology available in 5mm and 12mm Bladeless trocars, Dilating Tip trocars and Universal Sleeves. Also available in 12mm Blunt Tip trocars.



1. Comparing mean number of insertions for 5mm trocars out of 10 insertions with a quality index score of 92 or higher for 2B5LT (6.4 insertions) to COQ04/CFR03/CTRO3/CTF03/CFI13 (0.5 insertions), NB5STF/VS101005 (0.4 insertions) in benchtop tests, P<0.05. Comparing mean number of insertions for 12mm trocars out of 15 insertions with a quality index score of 92 or higher for 2B12LT (11.6 insertions) to CTF73/COR29/CFI83/CF773/CFR73/COR47 (0.4 insertions), NB12STF/VS101012P (1.6 insertions) in benchtop tests, P<0.05. Quality index score of 92 or higher indicates a less than 20% chance of having to rewipe the endoscope lens due to trocar induced smudging.

More options – more value

At Ethicon, we recognize that not all ports require the same level of trocar functionality. We offer the flexibility to choose the trocar for the right level of performance, when and where it's needed.

- ENDOPATH XCEL® Trocars with OPTIVIEW® Technology for camera ports
- ENDOPATH XCEL® Trocars for working ports
- ENDOPATH BASX® Trocars for retraction ports

Applicable specialties

- Bariatric • Colorectal • General • Gynecologic • Urologic

Bladeless Trocars

Code	Description	Length
2B5ST	Stability sleeve, 5mm diameter	75mm
2B5LT	Stability sleeve, 5mm diameter	100mm
2B5XT	Stability sleeve, 5mm diameter	150mm
2B12LT	Stability sleeve, 12mm diameter	100mm
2B12XT	Stability sleeve, 12mm diameter	150mm

Blunt Tip Trocars

Code	Description	Length
2H12LP	Stability sleeve, 12mm diameter	100mm

Instruments per Sales Unit: 6

Dilating Tip Trocars

Code	Description	Length
2D5ST	Stability sleeve, 5mm diameter	75mm
2D5LT	Stability sleeve, 5mm diameter	100mm
2D12LT	Stability sleeve, 12mm diameter	100mm
2D12XT	Stability sleeve, 12mm diameter	150mm

Universal Sleeves

Code	Description	Length
2CB5ST	Stability sleeve, 5mm diameter	75mm
2CB5LT	Stability sleeve, 5mm diameter	100mm
2CB12LT	Stability sleeve, 12mm diameter	100mm

How to Order

All purchase orders are made to Johnson & Johnson Health Care Systems, Inc. (JJHCS) or your authorized distributor.

Electronic Ordering Options

(Note: Placing orders electronically avoids minimum order fees for your hospital.)

Johnson & Johnson Gateway: www.jnjgateway.com/commerce

For questions about your order, visit the website or call 1-866-JNJ-GATE (565-4283).

Global Healthcare Exchange (GHX): www.ghx.com

For questions about your order, visit the website or call 1-800-YOUR-GHX (968-7449).

Electronic Data Interchange (EDI): JJHCS EDI Help Line: 1-800-262-2888

Nonelectronic/Manual Ordering Options

Call JJHCS at 1-800-255-2500 between 8:30 a.m. and 9:00 p.m. Eastern Standard Time or fax your order to 1-732-562-2212.

Customer Support

For product use assistance, clinical guidelines, service and repair, emergency assistance, copy of a 510(k) clearance letter or complaints, please contact our Customer Support Center at customersupport@eesus.jnj.com or by calling **1-877-ETHICON** (384-4266). Our Support Center is staffed 24 hours a day, 7 days a week by qualified nurses to answer your product-related questions.