


# ECHELON FLEX™ GST System

## Exceptional staple line integrity across the broadest range of tissue thickness<sup>1</sup>



<sup>1</sup> The ECHELON FLEX™ GST System was designed and tested to meet rigorous design requirements for staple line integrity. The GST System accommodates a tissue thickness range of 1.0mm to 4.0mm (measured at 8g/mm<sup>2</sup>, tissue comfortably compressed to closed staple height during firing per IFU) while the Medtronic Tri-Staple™ portfolio is intended for a tissue thickness range of 0.75mm to 3.0mm (per IFU & materials downloaded from Medtronic website on Nov 16, 2016).

The third-party trademarks used herein are trademarks of their respective owners.

©2017 Ethicon, Inc. All rights reserved. 081793-171005