


LIGACLIP®

Designed for **versatility**
and great **value**

LIGACLIP® Multi-Patient Use Single Clip Appliers



Features

- Lateral and transverse grooves on clip surface to help provide on-structure clip security
- Proximal to distal clip closure designed to provide minimal gap
- Grooves on the inside jaw surfaces of applicator designed to increase in-jaw clip security
- Tumbled tips to help reduce trauma to vessels
- Angled jaws designed to improve visibility beyond structures
- A mechanical stop in the LX100 & LX200 Series designed to limit the maximum amount of force applied to tissue.

ETHICON
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

Shaping
the future
of surgery

Product Overview

- LIGACLIP® Multi-Patient Use Single Clip Appliers provide a complete line of conventional, angled and right-angle products.
- The V-shaped clip is designed for ligation of tubular structures.
- The handles and clip cartridges are conveniently color-coded.
- Right-angle appliers include the LX130, LX230, LC330, LC430 models.

Ordering Information

Code	Description	Tip Angle Specification	Clip Size	# Clips	Clip Code
LX105	5-3/4", Small clips, Blue handle	15° angle	Small	6 20	LT100 ■ LT102 ■
LX107	7-1/2", Small clips, Blue handle	15° angle			
LX110	10-1/2", Small clips, Blue handle	15° angle			
LX130	10-1/2", right angled, Small clips, Blue handle	90° angle			
LX205	5-3/4", Medium clips, Silver handle	15° angle	Medium	6 20	LT200 ■ LT202 ■
LX207	7-1/2", Medium clips, Silver handle	15° angle			
LX210	10-1/2", Medium clips, Silver handle	15° angle			
LX220	10-1/2", angled jaw, Medium clips, Silver handle	45° angle			
LX230	10-1/2", right angled, Medium clips, Silver handle	90° angle			
LC307	7-1/2", Medium/Large clips, Green handle	15° angle	Med/Large	6	LT300 ■
LC310	10-1/2", Medium/Large clips, Green handle	15° angle			
LC320	10-1/2", angled jaw, Medium/Large clips, Green handle	15° angle			
LC407	7-1/2", Large clips, Yellow handle	15° angle	Large	6	LT400 ■
LC410	10-1/2", Large clips, Yellow handle	15° angle			
LC420	10-1/2", angled jaw, Large clips, Yellow handle	45° angle			
LC430	10-1/2", right angled, Large clips, Yellow handle	90° angle			

Instruments per Box: 1

How to Order

All purchase orders are made to Johnson & Johnson Health Care Systems, Inc. (JJHCS) or your authorized distributor.

Electronic Ordering Options

(Note: Placing orders electronically avoids minimum order fees for your hospital.)

Johnson & Johnson Gateway:
www.jnjgateway.com/commerce

For questions about your order, visit the website or call 1-866-JNJ-GATE (565-4283).

Global Healthcare Exchange (GHX):
www.ghx.com

For questions about your order, visit the website or call 1-800-YOUR-GHX (968-7449).

Electronic Data Interchange (EDI):
JJHCS EDI Help Line: 1-800-262-2888

Nonelectronic/Manual Ordering Options

Call JJHCS at 1-800-255-2500 between 8:30 a.m. and 9:00 p.m. Eastern Standard Time or fax your order to 1-732-562-2212.

Customer support

For product use assistance, clinical guidelines, service and repair, emergency assistance, copy of a 510(k) clearance letter or complaints, please contact our Customer Support Center at **customersupport@eesus.jnj.com** or by calling **1-877-ETHICON** (384-4266). Our Support Center is staffed 24 hours a day, 7 days a week by qualified nurses to answer your product-related questions.