

# CONTOUR<sup>®</sup> can. Reload after reload.

## CONTOUR<sup>®</sup> Curved Cutter Stapler



### Reloadable for multiple firings

- Designed to be reloaded when needed to complete a single transection with multiple firings
- Multiple firings may be used to prevent overloading of the tissue in the jaws of the device
- Knife extends to the retaining pin and staple line extends beyond the pin, supporting multiple firings
- Devices are reloadable and may be fired up to 6 times in a single procedure

### Designed for deep pelvis access

- Unique curved head designed to cut and staple deep in the pelvis
- Parallel closure and tissue retention pin designed to transect the rectum in a single perpendicular firing
- CONTOUR<sup>®</sup> allows placement of a 40mm cut line in the width of a 30mm space

**ETHICON**  
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

Shaping  
the future  
of surgery

# Product overview

- Four curvilinear rows of staples with a single cut between
- Allows low pelvic access without handle obstruction over existing linear staplers
- Tissue retaining pin with manual closure option holds tissue in place during compression and firing
- Parallel closure for controlled compression and less tissue movement
- Breakthrough washer provides audible and tactile feedback every firing



Device code	Tissue	Staple rows	Cut line length	Staple line length	Max firings	Qty/box
CS40B	■ Regular	4	40mm	51mm	6	3
CS40G	■ Thick	4	40mm	51mm	6	3



Reload code	Tissue	Closed staple height	Open staple height	Wire diameter	Qty/box
CR40B	■ Regular	1.5mm	3.5mm	0.20mm	6
CR40G	■ Thick	2.0mm	4.7mm	0.23mm	6

## Customer support

For product use assistance, clinical guidelines, service and repair, emergency assistance, copy of a 510(k) clearance letter, or complaints, please contact our Customer Support Center by calling **877-ETHICON (384-4266)**. Our support center is staffed 24 hours a day, 7 days a week by qualified nurses to answer your product-related questions.

# Ordering info

## How to order

All purchase orders are made to Johnson & Johnson Health Care Systems, Inc. (JJHCS).

If you want to order direct, you may order electronically (online) at:

- <https://us.jjcustomerconnect.com> or **1-866-565-4283**
- **Electronic Data Interchange (EDI)**— EDI Helpline: **1-800-262-2888**

Or, to place a non-electronic (manual) order, contact Johnson & Johnson Health Care Systems Inc. at 1-800-255-2500 between 8:30 a.m. - 6:30 p.m. (Eastern Standard Time) or fax us at **1-732-562-2212**.